## RITA® IntelliFlow™ Infusion Pump

**Caution:** Federal Law (USA) restricts this device to sale by or on the order of a physician. Read all instructions carefully before use.

### Indications For Use:

The RITA® IntelliFlow Infusion Pump when used with the RITA® Tubing Set and reusable occlusion bed, is designed for the local delivery of saline through the irrigation ports of the RITA® electrosurgical device during radiofrequency ablation procedures with the RITA® 1500x RF Generator with software version 8.10 or higher.

The RITA® IntelliFlow Pump is not to be used:

- for drug delivery
- for intravascular, intra-arterial or epidural delivery of fluids.
- during MRI procedures.

### Warning:

- Do not operate the IntelliFlow Pump in the presence of flammable anesthetics, oxygen-enriched or explosive atmospheres.
- If four or more tubes become occluded, improper or unpredictable lesion size may result

### Precautions:

- Do not use IntelliFlow Pump or RF Generator if either of them have been dropped or damaged.

### Equipment List:

- Sterile 0.9% normal injectable saline solution (at least 100 cc bag required).
- RITA® Tubing Set.
- IntelliFlow Pump
- Reusable Occlusion Bed
- Communication Cable, (700-102785)

### Front Panel Display

**"Power" LED** - A green light on the front of the panel indicates the power is on.

**"Status" LED** - A tri-color light-emitting diode (LED) on the front of the panel is used to indicate the pump drives condition. LED color versus status is as follows:

- Flashing Red - Call customer service
- Amber (Brownish orange) - Pump head tubing set not seated.
- Red - Communication error- Call customer service.
- Flashing Green - Waiting for commands from generator.
- Green - Ready for operation.

### Rear Panel

- "Power ON/OFF" Button. Switching this button turns the Pump power ON and OFF.
- "Power Cord" Connection. Port for connecting the RITA® Pump to power outlet via the power cord.
- "Serial Port" Connection. Port for connecting to Port B of the RF Generator. For use with StarBurst XLI or XLI-enhanced.*

### Instructions For Use:

- Connect the AC Power cord to pump then plug in to wall power.
- Using the Communication Cable, connect the IntelliFlow Pump to the RITA® RF Generator (Port B).
- For instructions on loading and purging the RITA® Tubing Set, please refer to the Instructions for Use packaged with the RITA® with Tubing Sets
- Turn pump on.
- Complete Purge process by pressing the "A" (Purge Mode) on the RF Generator. The Purge Mode will infuse fluid through system at a rate of 0.7 ml/min. for 300 sec.
- Purging is complete when there is continuous flow from all micro infusion needles.
- If the system is purged before the 300 seconds, press the "B" on the generator to stop the purging process.
- For pump operation and use of the entire system refer to the User's Guide and Service Manual for the RITA® 1500x RF Generator and the instructions for use for the RITA® Starburst XLI or XLI-enhanced.

IntelliFlow Pump with Tubing Set

### Maintenance

The IntelliFlow Pump is designed for indoor use in a dry operating room/procedure room environment. The IntelliFlow Pump requires **no maintenance or calibration** by the user. The user should not remove the cover. Removal of unit cover voids the warranty. All maintenance and calibration must be referred to Customer Service at AngioDynamics, Inc.

### Cleaning and Disinfecting the Pump

Shut off the IntelliFlow Pump, and disconnect power before cleaning. The IntelliFlow Pump should be given reasonable care and be kept clean and sanitary.

The IntelliFlow Pump may be cleaned with a damp disposable wipe moistened with 70/30 isopropyl alcohol and water. Do not use abrasives, caustics, or mineral spirits. Use of these agents to clean the IntelliFlow Pump, or any of its accessories may cause damage and voids the warranty. All electrical connections must be air-dried before use.

### Warning:

Electric shock hazard. Do not saturate the IntelliFlow Pump with liquids. Do not allow liquids to run inside the unit. Do not immerse the IntelliFlow Pump in water.

### Trouble Shooting Guide

The pump does not recognize generator

- Ensure pump is plugged into serial port B.

The pump does not turn on:

- Ensure power cord is plugged in and power turned on.

IntelliFlow Pump Specifications	
<b>Power Input</b>	100VAC-240VAC, 50/60hz
<b>Operating Temperature</b>	15°C - 35°C
<b>Storage Temperature</b>	-10°C - 65°C
<b>Humidity</b>	10% - 90%
<b>Flow Rate</b>	0.05 - 0.7 ml/min
<b>Max. Load</b>	60 psi
<b>Regulation (Line, Load, Drift)</b>	± 0.25%
<b>Dimensions</b>	14.75" x 8.8" x 6.25" (width x depth x height) (37.5 cm x 22.4 cm x 15.9 cm)
<b>Weight</b>	10 lbs. (4.5 kg)
<b>Controls</b>	Power on/off

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© 2009, AngioDynamics, Inc.
RITA e IntelliFlow son marcas comerciales de AngioDynamics, Inc.
RITA e IntelliFlow er varemærker, der tilhører AngioDynamics, Inc.
AngioDynamics UK Ltd.
RITA e IntelliFlow zijn handelsmerken van AngioDynamics, Inc.